

PRODUCT SPECIFICATIONS FOR CAST ROLLS

Purpose: Immobilize and protect fractures and other injuries.

Description: Cast rolls are knitted fiberglass tape coated with a water activated resin that hardens to weight bearing strength within 20 minutes to protect a patient's injury.

Use: A cast roll is removed from a moisture proof pouch and dipped into room temperature water and the excess water squeezed from the roll. The roll is then applied over stockinet and under cast padding to immobilize and protect the patient's injury and then the resin is given time to form a hard protective shell.

Product:

Item	Description	Quantity
2759002	2" X 4 Yard Cast Coated Roll	10/Box
2759003	3" X 4 Yard Cast Coated Roll	10/Box
2759004	4" X 4 Yard Cast Coated Roll	10/Box
2759005	5" X 4 Yard Cast Coated Roll	10/Box

Labeling:

Each cast roll is labeled with the size of the cast roll, color dot if applicable and the manufacturing date. Each box of 10 rolls is labeled with size, color dot(if other than white) and manufacturing date.

Packaging:

Each cast roll is packaged in a foil moisture proof pouch with ten pouches per box. Application instructions are included with each box.

Storage:

The cast rolls should be stored in a cool dry location.

Recommended storage temperature: 65F- 75F/ 18C- 24C